


Cy3 Azide

Product Number:	CCR-AZ119
Product Name:	Cy3 Azide
CAS Number:	N/A
Molecular Weight:	850.92
Solubility:	Water, DMSO, DMF, MeOH
Storage:	Upon receipt store at -20°C.
Appearance:	Red solid
Shelf Life:	At least 12 months at -20°C


Chemical Structure:



Optical Properties

λ_{abs}	550 nm
λ_{abs}	570 nm
ϵ_{max}	150,000 M ⁻¹ cm ⁻¹
Φ_f	0.10
CF ₂₆₀	0.05
CF ₂₈₀	0.05
Flow cytometry laser line	532 or 555 nm
Microscopy laser line	532 or 555 nm
Spectrally Similar Dyes	Alexa Fluor 555 CF Dye 555, DyLight 555

Absorption and Emission Spectra



For Research Use Only