


## Biotin-PEG4-NHS Ester

Product Number: CCR-B103

Product Name: Biotin-PEG4-NHS Ester

Chemical Structure:



CAS Number: 459426-22-3

Molecular Weight: 588.67

Chemical Composition:  $C_{25}H_{40}N_4O_{10}S$

Solubility: DMSO, DMF, DCM, Chloroform, THF

Storage: Upon receipt store at  $-20^{\circ}C$ . Product shipped at ambient temperature.

Appearance: White amorphous solid

Shelf Life: At least 12 months at  $-20^{\circ}C$

**For Research Use Only**