


## Alkyne-PEG4-NHS Ester

Product Number: CCR-TA103  
Product Name: Alkyne-PEG4-NHS Ester  
Alternative Name: Acetylene-PEG4-NHS Ester

Chemical Structure:



CAS Number: 1393330-40-9  
Molecular Weight: 401.41  
Chemical Composition: C<sub>18</sub>H<sub>27</sub>NO<sub>9</sub>  
Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.  
Appearance: Slightly yellow to orange oil at room temperature.  
Shelf Life: At least 12 months at -20°C

**For Research Use Only**