


## Biotin-PEG4-Alkyne

Product Number: CCR-TA105-25

Product Name: Biotin-PEG4-Alkne

Chemical Structure:



CAS Number: 1262681-31-1

Molecular Weight: 457.58

Chemical Composition: C<sub>21</sub>H<sub>35</sub>N<sub>3</sub>O<sub>6</sub>S

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: White solid

Shelf Life: At least 12 months at -20°C

**For Research Use Only**