


## Propargyl-NHS Ester

Product Number:	CCR-TA111
Product Name:	Propargyl-NHS Ester
Alternative Name:	3-Prop-2-ynyloxy-propanoic acid 2,5-dioxo-pyrrolidin-1-yl ester
Chemical Structure:	 <p>The chemical structure shows a propargyl group (a terminal alkyne) attached to a propyl chain via an ether linkage. This propyl chain is further connected to a propanoic acid moiety, which is esterified to the nitrogen atom of a 2,5-dioxo-pyrrolidine ring (NHS group).</p>
CAS Number:	1174157-65-3
Molecular Weight:	225.20
Chemical Composition:	C <sub>10</sub> H <sub>11</sub> NO <sub>5</sub>
Storage:	Upon receipt store at -20°C. Product shipped at ambient temperature.
Appearance:	Colourless solid
Shelf Life:	At least 12 months at -20°C

**For Research Use Only**