


Propargyl-Maleimide

Product Number: CCR-TA113

Product Name: Propargyl-Maleimide

Alternative Name: 1-(Prop-2-yn-1-yl)-1H-pyrrole-2,5-dione

Chemical Structure:



CAS Number: 209395-32-4

Molecular Weight: 130.12

Chemical Composition: C₇H₅NO₂

Storage: Upon receipt store at -20°C. Product shipped at ambient temperature.

Appearance: Slightly yellow to orange solid

Shelf Life: At least 12 months at -20°C

For Research Use Only